

Anaerobic Culture Bottles and Caps

Product Information

These culture bottles with 45 mm screw caps and 20 mm septum stoppers from Bellco Glass are especially suited to culture anaerobic bacteria and to maintain oxygen-free conditions during sample harvesting. The bottles come with a GL45 screw cap with silicone o-ring and a blue butyl rubber septum stopper, which can be repeatedly punctured with a needle to facilitate easy access to the culture. The bottles are available from 100 to 1000 ml and are manufactured from Borosilicate glass. The screw cap is made of Polypropylene. All components are autoclavable.



Anaerobic Culture Bottles and Caps

- Type 1, Class A Borosilicate Glass
- **GL45 PP screw cap**, including **silicone o-ring** and blue **butyl rubber septum stopper** for maintaining anaerobic conditions during sample harvesting

| Cat. No. | Capacity (ml) | Pcs./Cs. |
|-------------|---|----------|
| 2048-50100 | 100 | 12 |
| 2048-50250 | 250 | 12 |
| 2048-50500 | 500 | 6 |
| 2048-51000 | 1000 | 6 |
| 2048-50000 | Anaerobic cap with silicone o-ring and 1 septum stopper | 1 |
| 2048-11800A | Replacement septum stoppers | 100 |

Excerpt of our product range from Bellco Glass

- ▶ Anaerobic Culture Tubes (Balch Type and Hungate Type) and Bottles
- ▶ Aseptic Filling Bell
- ▶ Bioreactors, Vessels and Accessories
- ▶ Caps and Closures (KAP-UTS™, Silicone Sponge, Stainless Steel, Pouring Rings, Plugs)
- ▶ Cloning Cylinder
- ▶ Cover Slips (Round, Rectangular, Square, Photo-etched)
- ▶ Erlenmeyer Flasks
- ▶ Fernbach Flasks
- ▶ Fermenter
- ▶ Flasks/Bottles (Anaerobic, Biometer, DeLong®, Erlenmeyer-, Fernbach-, Nephelo- Roux-, Sample-, Shake-, Snap-Seal, Trypsinizing-, Volumetric-)
- ▶ Heating Baths
- ▶ Hybridization Ovens and Bottles
- ▶ Incubators (also Roll-In for Roller Apparatus)
- ▶ Insect Cell Culture Systems
- ▶ Magnetic Stirrers, Bell-ennium™, „Back-to Basics” (Single- and Multi Positions)
- ▶ Overhead Drives and Accessories
- ▶ Pipettes, Pipette Cans and Containers
- ▶ Pipette Plugger and Cotton Roving
- ▶ Rocker and Shaker
- ▶ Roller Apparatus
- ▶ Roller Bottles (Glass)
- ▶ Roller Drums
- ▶ Shake Flasks
- ▶ Shaker (Orbital-, Reciprocating-, Linear-)
- ▶ Slides (Plain/Frosted End, Counting Slides: Reich-, Somatic-, Fluorescence)
- ▶ Spinner Flasks (large selection of different types) and Accessories
- ▶ Sykes-Moore Chamber
- ▶ Tissue Homogenizers (Dounce, Potter-Elvehjem, Tenbroeck) and Sieves
- ▶ Trypsinizing Flasks
- ▶ Tubes (Anaerobic Balch-/Hungate-, Centrifuge-, Leighton-, with Screw Cap, Gas Fermentation)
- ▶ Water Baths (Shaking, with one or two separate pans)

Complete product list and individual product information on request