



6-well BioPress™ Culture Plates (BF-3000C)

Culture plate for use with the Flexcell® Compression system.

- > Each well contains an acrylic piston used with a stationary platen (supplied with the Compression system, not included with the plates) to deform 3D tissue samples or cell seeded constructs using the FX-5000™ Compression system.
- > Ring foam holders keep samples confined to central region of each well.
- > All plates come pre-sterilized.
- ➤ Available in cases of 10 or 40 plates.




Figure 33. BioPress™ culture plate with acrylic pistons adhered to the silicone elastomer well bottom, foam sample holders (far left), and stationary platens (left of plate)