


Culture Slips®

Surface treated slides for use with Streamer® and FlexFlow™.



- Available in two sizes:
75 × 25 × 1,0 mm (for Streamer® or FlexFlow™),
75 × 24 × 0,2 mm (for FlexFlow™ only).
- 75 × 25 × 1,0 mm Culture Slips® are rimmed with a 1,0 mm wide PTFE border to help limit cell culture growth to the portion of the slip exposed to fluid flow.
- Low autofluorescence.
- Matrix-treated to promote cell attachment.
- Covalently bonded surfaces: Amino, Collagen (Type I or IV), Laminin (YIGSR).
- Delivered in sterile twin packs for one time immediate use.
- Available in cases of 6 and 36 culture slips.

75 x 25 x 1,0 mm for Streamer® or FlexFlow™

Cat. No.	Description
CS-U	Culture Slips® – Untreated
CS-A	Culture Slips® – Amino
CS-C	Culture Slips® – Collagen Type I
CS-C/IV	Culture Slips® – Collagen Type IV
CS-L	Culture Slips® – Laminin

75 x 24 x 0,2 mm for FlexFlow™

Cat. No.	Description
FFCS-U	Culture Slips® – Untreated
FFCS-C	Culture Slips® – Collagen Type I