

## FlexLink® for Tension or Compression (FX-6000TFL or FX-5000CFL)

### Run multiple regimens at one time!

- Upgrade your existing FX-6000T or FX-5000C with a second controller that links with an FX-6000™ system computer allowing the user to simultaneously operate two different regimens on two different baseplates.  
*(May require an additional vacuum pump or compressor depending on the intensity of the regimens!)*
- Expand the application area of your FX-6000T system with the possibility to apply Compression tests by adding a Compression FlexLink® FX-5000CFL.  
*(Requires also a compressor!)*
- Expand the application area of your FX-5000C system with the possibility to apply Tension tests by adding a Tension FlexLink® FX-6000TFL.  
*(Requires also a vacuum pump!)*
- Connects in a series to operate up to four FlexLinks® from a central computer.



Figure 12. FX-6000T system with an additional FlexLink®