Junn

Osci-Flow[®] Flow Controller (STR-4000-OFS)

The ultimate in fluid flow direction control!

- Provides regulated oscillatory and pulsatile flow control via computer controlled action.
- > Minimizes flow response lag by eliminating inertial effects of decelerating and accelerating pumps and motors.
- Reverses fluid flow instantaneously.
- ➤ Works with Streamer[®] and FlexFlow[™] shear stress devices.
- Adapts to other perfusion systems.
- > Accomodates any MasterFlex L/S series or comparable tubing.
- Interfaces with most laptops via USB.
- Osci-Flow[®] System includes:
 Osci-Flow[®] device

 - Tubing and quick disconnects
 - StreamSoft[™] software



Figure 27. Osci-Flow[®] Controller




Figure 28. Set-up of the Streamer[®] device with Osci-Flow®