


## StageFlexer® Membranes

**Flexible growth surface for the StageFlexer® and FlexFlow™ devices.**

- Individually packaged 43 mm silicone rubber membranes (come in sterile culture dishes).
- Clear for direct viewing of cells optically.
- Viscoelastic under mechanical loading.
- Same growth surfaces as the BioFlex® culture plates.
- Covalently bonded surfaces: Amino, Collagen (Type I or IV), Laminin (YIGSR).
- Available in cases of 6 and 36 membranes.



Cat. No.	Description
SFM-U	StageFlexer® Membrane – Untreated
SFM-A	StageFlexer® Membrane – Amino
SFM-C	StageFlexer® Membrane – Collagen Type I
SFM-C/IV	StageFlexer® Membrane – Collagen Type IV
SFM-L	StageFlexer® Membrane – Laminin