

Flexcell® Transwell® Holder

Holds Transwells® to allow for cell migration and co-culture assays to be performed in Flexcell® culture plates.

- Analyze cell migration in response to strain.
- Available for 6- and 24-well Flexcell® culture plates.
- Compatible with standard Transwell® sizes for 6-well and 24-well culture plates.




Figure 9. 6-well Transwell® Holder.





Figure 10. Diagram of a Transwell® Holder with a Transwell® supported above the rubber membrane in the well of a BioFlex® plate

Transwell® is a registered trademark of Corning® Inc.