



NEU: Stapelbare Niedrigprofilplatten

Die neuen stapelbaren Niedrigprofil 96-Well Platten von Porvair Sciences sind durch ihre Flexibilität, Skalierbarkeit und dem großen Arbeitsvolumen perfekt für die Automation geeignet. Die Höhe der Platte von nur 27 mm reduziert den Platz gegenüber konventionellen 44 mm hohen Platten und macht diese ideal sowohl für manuelle Anwendungen als auch die Nutzung mit automatischen Systemen.

Spezifikationen:

- 27 mm hohe, stapelbare 96-Well Platte
- Bis zu 0,5 ml Arbeitsvolumen
- V-Boden und runde Wells
- Hohe Beständigkeit gegenüber Chemikalien und weiter Temperaturbereich von -196 °C bis +120 °C
- DNase- und RNase-frei
- Hergestellt aus Polypropylen
- Autoklavierbar und zentrifugierbar (4.800 x g)



- Alphanumerische Beschriftung der Wells
- Erhöhter Wellrand für einfaches Versiegeln z. B. mit den Sealern und Hitzeversiegelungsmatten von Porvair Sciences
→ Kontaktieren Sie uns für weitere Informationen!
- ANSI/SLAS Standard



Kat. Nr.	Beschreibung	VE
219007	96-Well stapelbare Niedrigprofilplatte, 27 mm, PP	80
219007FP	96-Well stapelbare Niedrigprofilplatte, 27 mm, PP	34 x 80
219017	96-Well stapelbare Niedrigprofilplatte, 27 mm, PP, steril	80
219017FP	96-Well stapelbare Niedrigprofilplatte, 27 mm, PP, steril	34 x 80

NEW: Low Profile, Stackable Microplates

The new stackable low profile 96-well plates from Porvair Sciences are ideal for automation due to flexibility, scalability and plenty of storage. By reducing the height of the plate from 44 mm conventional plates to 27 mm, the plates are ideal both for manual applications and for use with robotic systems.

Features:

- 27 mm high stackable 96-well plate
- Up to 0.5 ml working volume
- V-bottom and round wells
- High chemical and temperature resistance from -196 °C to +120 °C
- DNase- and RNase-free
- Manufactured from Polypropylen
- Autoclavable and possible to centrifuge (4,800 x g)



- Alphanumeric well coding
- Raised-wells for easy sealing e. g. with sealers and heat sealing foils from Porvair Sciences
→ Please contact us for further information!
- ANSI/SLAS standard

Cat. No.	Description	Cs.
219007	96-well stackable low profile plate, 27 mm, PP	80
219007FP	96-well stackable low profile plate, 27 mm, PP	34 x 80
219017	96-well stackable low profile plate, 27 mm, PP, sterile	80
219017FP	96-well stackable low profile plate, 27 mm, PP, sterile	34 x 80